

MRC - HighResolution


MR compatible HD video camera

up to 1,280 x 960 pixels

in-bore camera

no interferences

compact size



- Sensor Type: digital CMOS sensor, global shutter, color or b/w
Resolution: 1,280x960 pixel @43 Hz, 1,280x720 pixel @60 Hz, or ... (adjustable)
Interface: Gigabit Ethernet
- Dimensions: 83.5 mm x 37 mm; height: 21 mm + lens
Various lenses: standard M12 mount
Low weight: approx. 160 g

1995-2015
20 Jahre / 20 years
Innovation

MRC Systems GmbH - Hans-Bunte-Str. 8-10
D-69123 Heidelberg - Germany
info@mrc-systems.de - www.mrc-systems.de
phone: +49 (0) 6221 1380-300 - fax: -301

MRC